



MARSDEN
ACCURACY ASSURED

M-400

Marsden Portable Baby Scale

Lightweight and Portable

Ideal for Weighing Babies and Toddlers

Secure Hygienic Slide On/Off Weighing Pan

Ideal for Health Visitors and Midwives



MARSDEN
ACCURACY ASSURED



M-400

Description

Compact, stable and robust portable digital baby scale with secure hygienic slide on/off weighing pan, ideal for baby or stand-on toddler weighing.

Ideal for health visitors and midwives, this scale is easily transported, in either an easy to carry well padded bag or rucksack for 'hands free' operation.



Features

Capacity	20kg
Graduations	5g<10kg>10g
Dimensions	320mm x 310mm x 60mm
Weight of scale	3.3kg
Power Supply	6 x AA Batteries / Mains



Reliable. Accurate. Elegant.

