



MARSDEN
ACCURACY ASSURED

M-550

Marsden Low Cost Portable Floor Scale

Digital Floor Scale

Optional Bluetooth

Telehealth Applications

Ideal for GP Surgeries

Suitable for home use

Lightweight & portable



MARSDEN
ACCURACY ASSURED



M-550

Description

The Marsden M-550 Class III floor scale offers a real alternative for the mechanical scales in GP surgeries. It has a weighing accuracy of just 200g, which is five times better than the weighing performance of the mechanical scales. It has a wide weighing platform suitable to accommodate the larger patient.

The M-550 is low cost, but critically Class III compliant. It meets all the legal requirements of weighing patients, which makes it suitable to monitor the weight of patients in the surgery or at home. It is ideal for tele health programmes and has a Bluetooth option to connect the scales with other medical devices or databases.



CC-420

Features

Capacity	160kg
Graduations	200g<100kg>500g
Weight of scale	4kg
Power Supply	Mains / Batteries

Reliable. Accurate. Elegant.

