
Folding Privacy Screens



Folding Privacy Screens with Castors



Wall-Mounted Folding Privacy Screens



Privacy Screens

Shown: PV-3005



Shown: PV-3002



Folding Privacy Screens

The folding privacy screens are indispensable tools in every medical practice or facility.

They are available

- Without castors,
- With castors,
- Pull-out/wall-mounted.

Features

- Sturdy **profiled aluminium tubular frame** for enhanced **stability**.
- Individual panels joined by two **galvanized hinges** which can be folded in both directions.

- **Light-weight panels** made of 3 mm-thick, impact resistant and easy-to-clean **PVC rigid foam boards**.
- **Micro-cell foam structure** uniformly coloured, **flame-retardant** and resistant to repeated disinfectant use.

Options

The privacy screens panels are available in your **choice of the following colours at no extra charge:**

- White (similar to RAL 9003).
- Blue (similar to RAL 5017).
- Grey (similar to RAL 7037).
- Black (similar to RAL 9004).
- Yellow (similar to RAL 1018).
- Red (similar to RAL 3000).
- Green (similar to RAL 6024).

- **Mirror-like panels** also available (price **per panel** under *Spiegel-PV*).

- **All panels**, whether coloured or mirror-like, **can be combined and arranged individually** (please specify upon ordering).

Shown: PV-3004



Shown: PV-3003



Shown: Screen panel colours





Figure	Article Number	Panels	Dimensions (mm)		Weight (kg)	Price EUR
			Extended W x H	Folded W x H x D		
2	PV-3002	2	1140 x 1652	605 x 1652 x 56	7,70	
4	PV-3003	3	1720 x 1652	605 x 1652 x 82	11,60	
3	PV-3004	4	2260 x 1652	605 x 1652 x 105	15,50	
1	PV-3005	5	2840 x 1652	605 x 1652 x 120	19,40	
5, page 79	Spiegel-PV					

Privacy Screens

Shown: PV-3005/LR



Shown: PV-3002/LR



Folding Privacy Screens with Castors and Wall-Mounted Folding Privacy Screens

- The /LR folding privacy screens are equipped with castors, Ø 60 mm.
- The two outside panels are fitted with **rotating aluminium crossbars**. Each crossbar has one locking and one non-locking castor.
- The **inner panels** are also fitted at their junction with one **non-lockable castor**.

Options

The folding privacy screens with castors (**total height 1750 mm**) are also available for **wall-mounting** on one side.

The /WM wall-mounted folding privacy screens include:

- An aluminium rail with screws and dowels.


Shown: PV-3003/LRWM



Shown: PV-3003/LRWM



Shown: Hinge



Shown: PV-3003/LR, Spiegel-PV



Figure	Article Number	Panels	Dimensions (mm)		Weight (kg)	Price EUR
			Extended W x H	Folded W x H x D		
2	PV-3002/LR	2	1140 x 1750	605 x 1750 x 56	8,65	
5	PV-3003/LR	3	1720 x 1750	605 x 1750 x 82	12,65	
not shown	PV-3004/LR	4	2260 x 1750	605 x 1750 x 105	16,65	
1	PV-3005/LR	5	2840 x 1750	605 x 1750 x 120	20,65	
similar 3+4	PV-3002/LRWM	2½	1460 x 1750	605 x 1750 x 92	11,40	
3+4	PV-3003/LRWM	3½	2040 x 1750	605 x 1750 x 118	15,40	
similar 3+4	PV-3004/LRWM	4½	2580 x 1750	605 x 1750 x 131	19,40	
similar 3+4	PV-3005/LRWM	5½	3160 x 1750	605 x 1750 x 146	23,40	