



EMS Physio ... *making it better*

## EMS Physio ... *making it better*



James Greenham  
*Managing Director*  
[jgreenham@emsphysio.co.uk](mailto:jgreenham@emsphysio.co.uk)

Welcome to the 2015 edition of our product catalogue and thank you for choosing EMS Physio as a valued supplier.

Having recently celebrated our 90th year of trading, we continue to progress and improve our product range.

2015 sees the launch of a new and exciting generation of product software.

The user's experience with the Megapulse Senior shortwave unit is now greatly enhanced and improved with the new interface, with clear on screen instructions, colour touchscreens, informative pictures and helpful guidance notes.

Our core values continue - excellent equipment, supported by trained and knowledgeable staff, complemented by the highest quality and on-going service. Our global presence will reassure you that EMS Physio offers you the stability and support you need and expect from one of the world's leading physiotherapy brands.

As always, our product range continues to change, innovate and inspire you to offer the best technology for the effective treatment of your patients. Whether it's ultrasound, shortwave diathermy, couches, electrotherapy, laser, exercise equipment, gait analysis, shockwave therapy, compression therapy, EMG biofeedback or real-time medical imaging, EMS Physio remains your obvious choice.

Thank you for your commitment to EMS Physio.

James Greenham *Managing Director*



Cathy Camden *Export Manager*  
[+44 \(0\)1235 772 272](tel:+44201235772272) | [sales@emsphysio.co.uk](mailto:sales@emsphysio.co.uk)



Maya Eadie *Marketing Manager*  
[+44 \(0\)1235 772 272](tel:+44201235772272) | [sales@emsphysio.co.uk](mailto:sales@emsphysio.co.uk)



Adrian Keep *Area Sales Manager for South and S. Wales* | [07770 860 301](tel:07770860301) | [akeep@emsphysio.co.uk](mailto:akeep@emsphysio.co.uk)



Roy Bennett *Area Sales Manager for North, Midlands and N. Wales* | [07770 860 308](tel:07770860308) | [rbennett@emsphysio.co.uk](mailto:rbennett@emsphysio.co.uk)

## Pulse Press

Dynamic Compression Therapy has been clinically proven to significantly improve healing time in the treatment of oedema and other venous and lymphatic disorders of the lower limb. When used post surgery, reduces the chance of a D.V.T occurring.

### Physio 3



PCT030

### Physio 6



PCT031

### Physio 12




PCT033

### Lymph 12 Pro



PCT064

## Garment Sizes Key



### Technical Specification

**PHYSIO 3**  
**3 chamber sequential gradient compression**  
 Pressure 30-100 mmHg,  
 Compressor 14-20 ltrs/min  
 Cycle Times Inf. 5-120 secs, Def. 5-60 secs  
 Garments 1 ch and 3 ch  
 Size 166 x 135 x 150 mm  
 Weight 1.85 kg

**PHYSIO 6**  
**6 chamber sequential gradient compression**  
 Pressure 30-100 mmHg  
 Compressor 14-20 ltrs/min  
 Cycle Times 5-120 secs inflation, 5-60 secs deflation  
 Garments 3 ch and 6 ch  
 Size 187 x 155 x 150 mm  
 Weight 2.10 kg

**PHYSIO 12**  
**12 chamber sequential gradient compression**  
 Pressure 30-100 mmHg  
 Compressor 24-30 ltrs/min  
 Cycle Times 5-120 secs inflation, 5-60 secs deflation, 5-120 mind therapy  
 Garments 6 ch and 12 ch  
 Size 187 x 155 x 150 mm  
 Weight 3.15 kg

**LYMPH 12 PRO**  
**12 chamber sequential gradient compression**  
 Pressure 30-100 mmHg  
 Compressor 16-24 ltrs / min  
 Cycle Times Inflation Automatic, Deflation Automatic, Therapy 30-60 mins  
 Mode Selection x7 different modes  
 Size 187 x 166 x 150 mm  
 Weight 3.15 kg

### Ordering Information

Pulse Press Physio 3	PCT030
Pulse Press Physio 6	PCT031
Pulse Press Physio 12	PCT033
Lymph 12 Pro	PCT064

### Accessory Ordering Information

Garment to use with Physio 3	Garment Size Guide (cm)				
	A	B	C	D	E
Half Arm Sleeve	50	43	42		
Wide Long Arm Sleeve	60	76	47		
Half Legging	60	76	47		
Long Legging	70	74	32	31	50
Long Arm + Shoulder	60	76	47		
Long Leg	78	67	34	31	60
Ex. Wide Long Leg	90	75	40	31	80
Wide Long Leg	80	83	32	31	54
Wide Long Leg	84	90	32	31	53

Garment to use with Physio 6	Garment Size Guide (cm)				
	A	B	C	D	E
Wide Long Arm Sleeve	60	76	47		
Long Legging	70	74	32	31	50
Wide Long Legging	80	83	32	31	54
Wide Trouser Garment	130	80	120/13778/95		
Long Arm + Shoulder	60	76	47		
Ex. Wide Long Leg	90	75	40	31	80
Wide Long Leg	84	90	32	31	53

Garment to use with Physio 12 and Lymph 12 Pro	Garment Size Guide (cm)				
	A	B	C	D	E
Wide Trouser Garment	130	30	120/13778/95		
Wide Long Arm	60	76	47		
Long Arm + Shoulder	60	76	47		
Wide Long Leg	80	83	32	31	54
Bariatric Legging	92	82	40	31	80
Lymph Legging	82	80	120/13778/95		