

M-300

Marsden Portable Baby Scale

Ultra Lightweight

Portable, Ideal for Health Visitors
and Midwives

Optional Rucksack or Carry Case

Large, Clear Display

Suitable for Neonatal Baby Units

600 Hours Battery Life



M-300

Description

At 3kg it is genuinely lightweight & portable, ideal for health visitors & midwives. It can be supplied with a thoughtfully made rucksack or carry-case to enable easy carriage. It is one piece scale, so avoiding the need to add or remove the baby tray. It has a large, clear display. Its 2g accuracy makes the scale suitable for neo natal baby units within hospitals. The picostrain technology and new PCB design facilitate the accuracy and an extreme long battery life, at least 600 hours.



Features

Capacity	15kg
Graduations	2g<6kg>5g
Dimensions	370mm x 590mm
Weight of scale	3.1kg
Power Supply	2 AA Batteries / Optional Mains Adaptor

Reliable. Accurate. Elegant.

