



# EDAN

## PADECG A Mobile ECG Solution

PADECG consists of a professional ECG software and a wireless ECG acquisition module. It provides a mobile solution to meet the demands of caregivers in a mobile treatment environment. EDAN's reliable SEMIP algorithm is also provided in it, which allows you to have a fast and confident ECG assessment before sending the ECG report to Smart ECG Viewer for further analysis.




- A Mobile ECG Solution
- User-friendly Interfaces
- Extensive Communication



Care for Health



# PADECG Paperless Workflow



## ▪ Touch Screen

Operating PADECG is just as simple as using your smart phone. Convenient gestures such as long pressing, un-pinching, and dragging are all supported to enable users a better experience.

## ▪ ECG Sampling

Both pre-sampling and real-time sampling modes are offered to cater to different situation, no matter it is an emergent case, general ward rounds, or home visiting.

## ▪ Filter

Armed with the comprehensive filters, you can acquire superb ECG traces with full confidence just within several touches.

## ▪ Algorithm

With EDAN's professional ECG algorithm SEMIP, PADECG gives accurate diagnosis results and offers doctors a reliable reference.



### Minimum Tablet Configuration requirement

PADECG					
Screen size	7"	Battery	4,000 mAh	OS	Android 4.0 OS
Resolution	1024*600 (WSVGA)	Chipset	1.2GHz Dual Core	Memory	16GB + 1GB RAM
Camera	With Auto-focus function	Wi-Fi	Yes	Bluetooth	Yes

**EDAN** Edan Instruments, Inc.

- 3/F – B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
- Tel: +86-755-26898326 ▪ Fax: +86-755-26898330 ▪ www.edan.com.cn ▪ Email: info@edan.com.cn
- All rights reserved. Features and specifications are subject to change without notice



CE 0123

ENG-ECG-PADECG-V1.0-20140325