

The Art of Medical Design



We know the art of medical design

Our goal at Fimet is to make your work easier. We work together with dental professionals from all over the world and bring those ideas and suggestions into our designs. We are a partner that reacts to your needs quickly and flexibly.

We are a Finnish based company but through our global distributor partners, we operate in more than 50 countries. Our development, design and manufacturing are based in Finland.

We use only the highest quality materials to ensure state-of-the-art design and experience. You can also give the Fimet dental system a personal touch with a wide selection of colors and patterns to choose from.

Sincerely, Jan J. Tollet President and CEO, Fimet Group 600



Custom Made Dental Treatment Systems for Dentists and Oral Health Care Specialists

Improving Your Patient Care Experience with a Smile

Why dental professionals choose FIMET over other brands

- Features and functions custom made to suit your requirements
- Quality made in Finland
- Easy to start using and affordable to maintain
- Provides maximum comfort for the patients during treatment
- Happy healthy patients repeated visits
- Helps dental professionals improve clinical outcomes and enhance productivity
- Cares for your clinic image

F1 Traditional

The Fimet F1 dental treatment system can be equipped also with hanging instrument hoses for those users who prefer this way of workflow.

4



Simply elegant

F1's elegant and stylish design includes many small details that have a big impact on the way you work.

Safe choice

Keep your patients safe. Smooth surfaces and seamless upholsteries make both cleaning and disinfection easy.

Designed to work

Adjust the chair position and instrument speed ranges so that the chair fits your way of working.

Compact

Small base means more working space for both the dentist and assistant. Everything is within easy reach.

F1 Continental

The Fimet F1 dental treatment system combines innovative solutions with elegant design. The treatment systems are ergonomic, high-quality units that can be customized for each user's needs and preferences.

Simply elegant

F1's elegant and stylish design includes many small details that have a big impact on the way you work.

Safe choice

Keep your patients safe. Smooth surfaces and seamless upholsteries make both cleaning and disinfection easy.

Designed to work

Adjust the chair position and instrument speed ranges so that the chair fits your way of working.

Compact

Small base means more working space for both the dentist and assistant. Everything is within easy reach.



For You and Your patient

The Fimet F1's modularity eases the workflow for the dentist and assistant. It also makes the treatment more pleasent for the patient. Control the system with a wireless foot control, joystick or user panel.



Height can be adjusted between 40 to 91 cm



Small base that fits even the smallest spaces



Adjustable chair position and instrument speed



Supports camera and display integration



Easy to use multifunctional wireless foot control



Arm rests rotate 90° sideways and are removable



Space for 5 instruments on the bridge

First-class features



Only well-known quality instruments are used in the F1 treatment systems.

The flexible suction arm is equipped with high and low volume suction tubings plus further options.



F1 has an integrated flushing system for all instruments using water.



The water unit can be equipped either with valves for wet suction system or a separator intended for dry suction system.





The arm rests can be turned away and even removed to allow easy access.



The wide movement range of the F1 patient chair makes working even in upright position possible.



F1 Cart

The Fimet F1 Cart can be easily positioned for the requirements of various operations, or even moved to another room. It can be used as a standalone unit or together with a F1 patient chair, suction arm and a cuspidor.



Just plug compressed air and power to the unit and it is ready for work.



The Art of Medical Design



F1 Ponychair saddle stool



F1 Dental stool

Stools

The F1 stools can be delivered in the same upholstery colors as the F1 patient chair.





F1 Classic stool

Colors

Seamless, easy-to-clean and comfortable upholsteries are available in many fashionable colors.

F1 upholstery colors
• synthetic leather, PVC

F1 Royal Blue F11 Black F19 F21 Midnight Blue Anthracite Griffine • synthetic leather, PVC NG02 NG01 Black Blue NG04 NG03 Purple Violet

The specifications and colors in this leaflet may vary from the actual product and are subject to change without prior notice.



Options and features

Fimet F1 treatment system can be configured with a large variety of properties and options according to your wishes and preferences.

F1 Instrument Bridge
Space for five instruments and can include the following options:
Electric micromotor(s) with or without brushes, with light
Turbine tubing(s) with light
Ultrasonic scaler, with or without light
3- or 6-way syringe
Polymerization light
Intraoral camera

Suction Unit		
Includes:		
High and low volume suction	oses	
Optional syringe		
Optional polymerization light		
Optional intraoral camera		

Water Unit
Can be equipped with:
Integrated amalgam separator
Valves for wet suction system
Water bottle
Detachable porcelain bowl
Suction filter
Water heater
-

Standards
CE-approved according to Medical Device Directive 93/42/EC
IEC 60601-1: 2005 tested
ISO 9001 and ISO 13485 certified quality system

F1 Patient Chair Includes: Arm rests rotate 90 degrees to side to allow easy access, also removable Four individually adjustable settings including last-position-program Two-joint adjustable head rest Knee-break chair

Further features and options:

Led patient light
Support for intraoral camera and display integration
Wireless, multifunctional foot control for chair and instrument functions
Tray holder for one or two standard trays



Our mission is to design and produce affordable, high-quality dental treatment systems to meet your needs and expectations.

15

FIMET

Fimet Oy Teollisuustie 6 FI-07230 Monninkylä Finland

Tel: +358 19 521 6600 Fax: +358 19 521 6666 fimet@fimet.fi www.fimet.fi